



## Short report

## Study of DXS9895 and DXS7130: Population data from North of Portugal

Raquel Carvalho MSc, PhD Student<sup>a,\*</sup>, M. Fátima Pinheiro PhD, Forensic Expert<sup>b,c,d</sup><sup>a</sup> Oporto University, Porto, Portugal<sup>b</sup> National Institute of Legal medicine, North Branch, Portugal<sup>c</sup> Fernando Pessoa University, Porto, Portugal<sup>d</sup> ICBAS, Abel Salazar Institute for the Biomedical Sciences, Porto, Portugal

## ARTICLE INFO

## Article history:

Received 29 September 2009

Received in revised form

25 October 2010

Accepted 24 November 2010

Available online 8 December 2010

## Keywords:

North of Portugal

X-STRs

Population data

DXS9895

DXS7130

## ABSTRACT

The study of X-chromosomal short tandem repeats (X-STRs) can complement data obtained with autosomal and Y-STRs. This population study only concerns two X-STRs in order to add complementary data obtained with other X-STRs already studied by our laboratory. DXS9895 and DXS7130 were used to study a population sample of North of Portugal (101 female and 118 male samples). DNA was amplified in a multiplex reaction mix and the automatic detection was performed using capillary electrophoresis. Allele frequencies and several forensic parameters were calculated.

© 2010 Elsevier Ltd and Faculty of Forensic and Legal Medicine. All rights reserved.

## 1. Population

A total of 219 healthy and unrelated individuals (101 female and 118 male samples), from North of Portugal.

## 2. Methods

## 2.1. DNA extraction and amplification

DNA was extracted from blood samples using the Chelex<sup>®</sup> method<sup>1</sup> and purified, if necessary (when the sample failed to amplify in PCR reaction), using a modified organic phenol-chloroform-isoamylalcohol method. DXS9895 and DXS7130 were amplified in one multiplex with previously designed primers<sup>2,3</sup> using Qiagen Multiplex kit (Qiagen). PCR was performed in a total volume of 10 µL containing 5 µL of 2× Qiagen Multiplex PCR Master Mix (Qiagen), 1 µL of 10× Primer Mix, 3.5 µL of distilled water and 0.5 µL of template DNA. The final primer concentration in reaction mix was: DXS9895 (0.5 µM) and DXS7130 (4 µM). Thermocycling

conditions, using GeneAmp<sup>®</sup> PCR System 9700 (Applied Biosystems) were: an initial denaturation for 15 min at 95 °C, followed by 10 cycles of 30 sec at 94 °C, 90 sec at 60 °C, 60 sec at 72 °C and 20 cycles of 30 sec at 94 °C, 90 sec at 58 °C, 60 sec at 72 °C and a final extension of 60 min at 60 °C.

Aliquots containing 1 µL of amplified DNA were mixed with 13.5 µL of Hi-Di Formamide (Applied Biosystems) and 0.5 µL of LIZ 500 (Applied Biosystems) as internal standard. After denaturation (5 min at 95 °C), samples were injected in ABI PRISM<sup>®</sup> 3100 Genetic Analyser (Applied Biosystems). Fragment sizes were determined automatically using GeneScan<sup>®</sup> Analysis Software v3.7 (Applied Biosystems) and genotyping was performed through comparison with DNA control samples 9948 (Promega), K562 (Promega) and 9947A (Applied Biosystems) according to the recommendations of Szibor.<sup>4</sup>

## 2.2. Data analysis

Allele frequencies, Hardy–Weinberg equilibrium in female samples (exact test) and linkage disequilibrium in male samples were calculated using GENEPOP version 3.4 software package<sup>5</sup> (Table 1). Several parameters of forensic interest were estimated with the formulae proposed by Desmarais<sup>6</sup> (Table 2). Genetic distance estimations based on the number of different alleles ( $F_{ST}$ ) were calculated using software ARLEQUIN v3.11.<sup>7</sup>

\* Corresponding author. Delegação do Norte do Instituto Nacional de medicina Legal, I.P., Jardim Carrilho Videira, 4050 – 167 Porto, Portugal. Tel.: +351 22 207 38 50; fax: +351 22 332 59 31.

E-mail address: [maria.raquel.carvalho@gmail.com](mailto:maria.raquel.carvalho@gmail.com) (R. Carvalho).

**Table 1**

Allele frequencies of DXS9895 and DXS7130 in a population sample of North of Portugal ( $n = 218$ ).

	DXS9895	DXS7130
10	—	0.003
11	—	0.038
12	0.003	0.111
13	0.251	0.035
13.3	—	0.035
14	0.222	0.006
14.3	—	0.190
15	0.346	—
15.2	0.003	—
15.3	—	0.410
16	0.146	—
16.2	0.010	—
16.3	—	0.149
17	0.013	—
17.2	0.006	—
17.3	—	0.016
18.3	—	0.006
p	0.706	0.137

p: Hardy–Weinberg equilibrium exact test in the female sample.

### 3. Results and other remarks

Allele frequencies are shown in Table 1. The population under study was in Hardy–Weinberg equilibrium and the two X-STR (DXS9895 and DXS7130) were not in linkage disequilibrium ( $p = 0.838$ ). Several forensic parameters are shown in Table 2.

Population differentiation between this population sample (North of Portugal) and Spanish population<sup>8</sup> was evaluated by genetic distance analysis. It was observed that there was no significant genetic distance between these two populations for the two loci under study (DXS9895,  $F_{ST} = -0.002$ ; and DXS7130,  $F_{ST} = 0.001$ ).

This population study concerns only two X-STRs (DXS9895 and DXS7130). However, these two markers used with other X-STRs

**Table 2**

Forensic parameters of DXS9895 and DXS7130.

	DXS9895	DXS7130
PIC	0.748	0.759
PD female	0.890	0.906
PD male	0.742	0.750
PE trio	0.703	0.728
PE motherless	0.564	0.595

PIC – polymorphism information content; PD – power of discrimination; PE – power of exclusion.

could be useful in forensic practice, particularly in “deficient paternity” and other kinship cases.

#### Conflict of interest

None.

#### References

- Walsh PS, Metzger DA, Higuchi R. CHELEX® 100 as a medium for simple extraction of DNA for PCR-based typing from forensic material. *Biotechniques* 1991; **10**:506–13.
- Edelmann J, Deichsel D, Hering S, Plate I, Szibor R. Sequence variation and allele nomenclature for the X-linked STRs DXS9895, DXS8378, DXS7132, DXS6800, DXS7133, GATA172D05, DXS7423 and DXS8377. *Forensic Science International* 2002; **129**:99–103.
- Edelmann J, Szibor R. The X-linked STRs DXS7130 and DXS6803. *Forensic Science International* 2003; **136**:73–5.
- Szibor R, Edelmann J, Hering S. Cell line DNA typing in forensic genetics: the necessity of reliable standards. *Forensic Science International* 2003; **138**:37–43.
- Raymond M, Rousset F. GENEPOP (version 1.2): population genetics software for exact tests and ecumenicism. *Journal of Heredity* 1995; **86**:248–9.
- Desmarais D, Zhong Y, Chakraborty R, Perreault C, Busque L. Development of a highly polymorphic STR marker for identity testing purposes at the human androgen receptor gene (HUMARA). *Journal of Forensic Science* 1998; **43**:1046–9.
- Excoffier L, Lischer HEL. Arlequin suite ver 3.5: a new series of programs to perform population genetic analyses under Linux and Windows. *Molecular Ecology Resources* 2010; **10**(3):564–7.
- Zarrabeitia MT, Alonso A, Martin J, Gonzalez-Gay MA, Martin-Escudero JC, Pancorbo MM, et al. Study of six X-linked tetranucleotide microsatellites: population data from five Spanish regions. *International Journal of Legal Medicine* 2006; **120**:147–50.